


WEBINAR

Practical Small Animal Ultrasound: Scanning the Pancreas and Adrenals for Common Pathologies

November 2021



Your Host



Dr. Oron Frenkel, MD, MS

Emergency Physician & POCUS Educator

Chairman, Clarius Medical Advisory Board


GI, Pancreatic and Adrenal Ultrasound

“This is another **new area** where ultrasound examinations delineate **morphological abnormalities**. With an improved resolution of recent ultrasound machines, the pancreas is **readily visible in normal dogs**.”


Takayoshi Miyabayashi. GI, Pancreatic and Adrenal Ultrasound. World Small Animal Veterinary Association World Congress Proceedings, 2001. Source: <https://www.vin.com/apputil/content/defaultadv1.aspx?pId=8708&id=3843781>



“Pancreatic diseases and abnormalities frequently investigated by means of ultrasonography include pancreatitis, pseudocysts, abscesses, neoplastic lesions, and nodular hyperplasia.”



“...ultrasonography is useful in diagnosing pancreatic disease, guiding aspirates and biopsies, and monitoring response to treatment.



Ultrasonography of the Adrenal Glands in the Dog, Cat, and Ferret


“Ultrasonography has become an **important modality** for imaging adrenal glands in small animals. Normal adrenal glands and structural abnormalities, such as changes in **size, shape, and echogenicity**, can be identified.”

Barthez PY, Nyland TG, Feldman EC. Ultrasonography of the adrenal glands in the dog, cat, and ferret. *Vet Clin North Am Small Anim Pract.* 1998 Jul;28(4):869-85. doi: 10.1016/s0195-5616(98)50082-4. PMID: 9698619. Source: <https://pubmed.ncbi.nlm.nih.gov/9698619/>



Feline Abdominal Ultrasonography: What's Normal? What's Abnormal?


“Ultrasound examination of the adrenal glands can provide **important information** pertaining to several conditions including **hyperaldosteronism** and **hyperadrenocorticism**.”



Griffin S. Feline abdominal ultrasonography: what's normal? what's abnormal? The adrenal glands. J Feline Med Surg. 2021 Jan;23(1):33-49. doi: 10.1177/1098612X20979509. PMID: 33403910. Source: <https://pubmed.ncbi.nlm.nih.gov/33403910/>

Poll

What do you see as the key advantages of using ultrasound for veterinary medicine?



Your Expert Guest Speaker



Dr. Camilla Edwards
DVM, CertAVP, MRCVS

*Peripatetic Veterinary Ultrasonographer |
Educator | First Opinion Veterinary Ultrasound*

Pancreas and the Adrenals

Dr Camilla Edwards DVM CertAVP MRCVS
First Opinion Veterinary Ultrasound

Conflict of interest declaration: Honorarium

Flat commission for ultrasound machines sold through my website

What will we learn in this webinar?

- How to find the adrenals and the pancreas
- What are the landmarks to help you consistently find these structures
- Indications for scanning these organs.
- Signs of pathology to watch out for in the adrenals and pancreas



What knowledge are we already assuming?

We're assuming you have some basic knowledge about ultrasound machine set up and basic anatomy

Adrenals




Adrenals indications for scanning

- Suspect or confirmed hyperadrenocorticism to differentiate between pituitary dependent and adrenal dependent disease
- Suspect or confirmed hypoadrenocorticism
- Part of the normal full abdominal ultrasound examination

What am I looking for?

- Cats
 - Length 0.89-1.25cm
 - Cranial pole thickness 0.3-0.48cm
 - Caudal pole thickness 0.3-0.45cm
- Dogs
 - Length 1-5cm
 - Caudal pole
 - <10kg <0.54cm
 - 11-30kg <0.68cm
 - >30kg <0.8cm
- Left - peanut shaped
- Right – L shaped/arrow head
- Phrenicoabdominal vein



Adrenals – ultrasound artifacts


Normal artefacts

- Refraction from curved surface of kidneys.
- Gas from colon and duodenum (particularly annoying on right adrenal!)

Signs of pathology

- Acoustic shadowing from mineralisation (abnormal in dogs, incidental in cats)


Left adrenal – how to scan





Adrenals – How to scan video left adrenal




Adrenals - anatomy



Left Adrenal



Left adrenal – normal ultrasound images



Left adrenal– right lateral recumbency

Follow Aorta cranial from trifurcation in longitudinal view (probe in cranial – caudal orientation)


A hooked artery will leave the Aorta at the level of the kidney – renal artery

If you see 2 anechoic circles these are the cranial mesenteric artery and the coeliac artery


Between the cranial mesenteric artery and the renal artery (often closer to the latter) fan.

Tip – you won't necessarily have all the landmarks in one image – but when you are close to these you will find the left adrenal



Right adrenal – how to scan





Adrenals – How to scan video right adrenal



Right Adrenal



Right Adrenal – normal images



Adrenal glands - common pitfalls

- Mistaking a blood vessel for an adrenal gland (use Colour Doppler if in doubt)
- Mistaking the left adrenal gland for the right adrenal gland (move further cranial)


Right adrenal – left lateral recumbency

Find the CVC by scanning dorsally pointing towards spine (probe oriented cranial-caudal) and fanning.

Follow the CVC cranially – you will see it disappear ventrally away from the Aorta which continues cranial in a dorsal position.

Where the CVC descends ventrally fan between Kidney, Aorta and CVC.

The right adrenal is often very close along the CVC and is arrow shaped.



Adrenals – Pathology

Hyperadrenocorticism

- ACTH Dependent (PDH) – 80% of cases – bilaterally normal or enlarged adrenal glands
- ACTH Independent – 20% of cases (adenoma, adenocarcinoma, phaeochromocytoma) – unilateral adrenomegaly and contralateral adrenal atrophy

Hypoadrenocorticism

- Smaller than normal

Focal enlargement



- Benign adenoma, adenocarcinomas, phaeochromocytomas
- Local vessel invasion (phrenicoabdominal vein and CVC)

Mineralisation


- Cats normal, dogs sign of malignancy

Left adrenal mass

- Rotweiller
- 6yo
- Left adrenal caudal pole 5cm invading local vasculature



Pancreas





Pancreas - indications for scanning

- Cranial abdominal pain
- Vomiting
- Anorexia
- Lethargy
- **Pancreatitis!**

Pancreas – what is normal?

- Size, Shape and Postition
 - Left limb (easier in cats)
 - Right limb (easier in dogs)
 - Body (hard in all!)
- Margins
 - Smooth margins
- Echogenicity
 - Isoechoic
 - Homogenous
 - Hypoechoic centre – pancreatic duct (left limb) & pancreaticoduodenal vein (right limb)
- Distribution of abnormalities
 - Try to view in two planes.
 - Important to register the mesentery/fat surrounding the pancreas for abnormalities


Pancreas – ultrasound artifacts

- Gas in the gastrointestinal tract –
duodenum, transverse colon, stomach

anatomy

This diagram illustrates the female reproductive system. The central organ is the pear-shaped uterus, colored pink. A fallopian tube extends from the lower side of the uterus, also in pink. To the right of the fallopian tube is a large, dark red, bean-shaped ovary. Several red and blue lines represent blood vessels connecting the ovary to the uterus. Below these organs is a layer of blue, wavy, scalloped tissue representing the endometrium or vaginal lining. The entire diagram is set against a black background.

Right Limb



Pancreas – How to scan right limb video



Right Limb Pancreas – normal ultrasound images



Right limb and body of the pancreas – left lateral recumbency

Find the caudal end of the right kidney


Pointing probe down towards table, slide ventrally

The most dorsal, superficial, straight loop of small intestine is the duodenum

The right limb of the pancreas lies close to the duodenum (medially)

Follow the duodenum cranially to the pyloroduodenal junction to see the body of the pancreas between the stomach and the portal vein

Left Limb



Pancreas – How to scan left limb video



Left limb Pancreas – normal ultrasound images



Left limb of the pancreas – right lateral recumbency

Find the colon dorsal to the bladder

Follow the colon cranially in longitudinal view until the transverse colon is reached


Aim to view the spleen, stomach and transverse colon in the same image

Fan in this location particularly following the splenic vessels down to the portal vein (anechoic tubes)

Pancreas pathology

- Pancreatitis
 - Hypoechoic pancreas
 - Hyperechoic peripancreatic structures (steatitis)
 - Free fluid
 - Increased thickness
 - Irregular margins

THE PANCREAS BECOMES MORE OBVIOUS WITH PATHOLOGY



How good is ultrasound at diagnosing pancreatitis?

J Vet Intern Med 2013;27:913–918

Ultrasonographic Findings of the Pancreas in Cats with Elevated Serum Pancreatic Lipase Immunoreactivity

J.M. Williams, D.L. Panciera, M.M. Larson, and S.R. Werre

Background: Pancreatitis is a common disease in cats that is difficult to diagnose.

Hypothesis/Objectives: To determine the sensitivity and specificity of ultrasonographic changes of the pancreas with serum feline pancreatic lipase immunoreactivity (fPLI) as the standard for diagnosis of pancreatitis.

Animals: 35 cats with clinical signs consistent with pancreatitis with an abdominal ultrasound examination and serum fPLI concentration measured within 3 days of the ultrasound.

Methods: Retrospective study: Pancreatic thickness, pancreatic margination, pancreatic echogenicity, and peripancreatic fat echogenicity were evaluated. Sensitivity and specificity were calculated with an elevated serum fPLI concentration indicative of pancreatitis as the standard for diagnosis.

Results: Serum fPLI was elevated and diagnostic for pancreatitis in 19 of 35 cats. The single ultrasound characteristic with the highest sensitivity was hyperechoic peripancreatic fat at 68% (95% confidence interval = 44–87%), indicating a moderate probability that cats with pancreatitis will have this abnormality on ultrasonographic examination. Specificity was >90% for each of increased pancreatic thickness, abnormal pancreatic margin, and hyperechoic peripancreatic fat. The sensitivity and specificity of ultrasound were 84% (95% confidence interval = 60–97%) and 75% (95% confidence interval = 48–93%), respectively, in cats with elevated serum fPLI indicative of pancreatitis.

Conclusions and Clinical Importance: The presence of a thick left limb of the pancreas, severely irregular pancreatic margins, and hyperechoic peripancreatic fat in cats with appropriate clinical signs and elevated serum fPLI are highly supportive of pancreatitis.

Key words: Abdominal ultrasound; Pancreatitis; Serum feline pancreatic lipase immunoreactivity.

Received: 31 August 2019 | Accepted: 19 December 2019

DOI: 10.1111/jvim.15693

STANDARD ARTICLE

Journal of Veterinary Internal Medicine 
American College of Veterinary Internal Medicine

Association between abdominal ultrasound findings, the specific canine pancreatic lipase assay, clinical severity indices, and clinical diagnosis in dogs with pancreatitis

Harry Cridge¹ | Alyssa M. Sullivant¹ | Robert W. Wills² | Alison M. Lee¹

¹Department of Clinical Sciences, College of Veterinary Medicine, Mississippi State University, Mississippi State, Mississippi

²Department of Basic Sciences, College of Veterinary Medicine, Mississippi State University, Mississippi State, Mississippi

Correspondence

Harry Cridge, Department of Clinical Sciences,
College of Veterinary Medicine, Mississippi
State University, Mississippi State, MS 39762.
Email: harry.cridge@gmail.com

Abstract

Background: A clinical diagnosis (CDx) of pancreatitis includes evaluation of clinical signs, abdominal ultrasound (AUS), and pancreatic lipase. However, practitioners are using AUS to diagnose pancreatitis and are using AUS severity to guide decisions. The validity of this is unknown.

Objectives: To determine whether (1) there is a correlation between AUS, specific canine pancreatic lipase (Spec cPL) assay, and CDx; (2) individual AUS abnormalities correlate more closely with CDx than others; (3) AUS severity mirrors clinical severity indices; (4) changes in AUS can be used as a marker for changes in Spec cPL or CDx; and (5) the sensitivity and specificity of AUS for pancreatitis.

Animals: One hundred fifty-seven dogs.

Methods: In this retrospective case study, inclusion criteria were signs of gastrointestinal, pancreatic disease, or both, in addition to having a Spec cPL and AUS performed within 30 hours. Information extracted from the records included bloodwork, Spec cPL, AUS images/clips, and severity of ultrasonographic findings.

Results: AUS was weakly correlated with Spec cPL ($r_s = .0178$, $P = .03$) and moderately correlated with CDx ($r_s = .379$, $P = <.001$). Pancreatic size ($r_s = .285$, $P = <.001$), echogenicity ($r_s = .365$, $P = <.001$), and mesenteric echogenicity ($r_s = .343$, $P = <.001$) were correlated with CDx. Change in AUS was not correlated with Spec cPL or CDx changes. When pancreatic enlargement, echogenicity, or altered mesenteric echogenicity were required for a diagnosis, the sensitivity and specificity were 89% (95% confidence interval [CI] 71.8, 97.7) and 43% (95% CI 34.0, 51.6). When all 3 criteria were required, the sensitivity and specificity were 43% (95% CI 24.5, 62.8) and 92% (95% CI 85.3, 95.7).

Conclusions: AUS should not be used in isolation to diagnose pancreatitis and is a poor indicator of severity.

Abbreviations: APACHE, acute patient physiologic and laboratory evaluation; AUS, abdominal ultrasound; CAPS, canine acute pancreatitis severity; CDx, clinical diagnosis; CI, confidence interval; ICA, ionized calcium; Spec cPL, specific canine pancreatic lipase; UPASS, ultrasonographic pancreatic assessment severity score.

This is an open access article under the terms of the Creative Commons Attribution-NonCommercial License, which permits use, distribution and reproduction in any medium, provided the original work is properly cited and is not used for commercial purposes.
© 2020 The Authors. *Journal of Veterinary Internal Medicine* published by Wiley Periodicals, Inc. on behalf of the American College of Veterinary Internal Medicine.

What are the take home messages today?

- **Know your landmarks**
 - **RLP** – duodenum
 - **LLP** – stomach, transverse colon, spleen
 - **Body** – Stomach
 - **Left adrenal** – Aorta, renal artery, cranial mesenteric artery, coeliac artery
 - **Right adrenal** – Caudal vena cava and Aorta divergence
- **What changes are you looking for?**
 - Pancreas – hyperechoic fat, size, margins, free fluid.
 - Adrenal glands - change in size in caudal pole thickness
 - Compare the size to the contralateral gland.
 - Vessel invasion
- **PRACTICE FINDING THE LANDMARKS!**



Thank you!

- Get in touch with any questions
Camilla@fovuc.co.uk
- Check out my website
www.fovuc.co.uk
 - Scanning
 - Reviews
 - Courses
 - Community

Live Demonstration

5 Ways POCUS Can Help in an Emergency

1. *Superior to radiography for the detection and scoring of free intraabdominal fluid.*
2. *Real-time information for rapid diagnosis, prompting life-saving maneuvers.*
3. *Cageside and low impact with minimal restraint.*
4. *Safe, no shaving required and radiation sparing.*
5. *App-based POCUS is easier to learn and use.*



Adam Behrens, VMD

Wandering Vet

“ The Clarius C7 HD microconvex scanner provides images as good as I’ve seen on large console machines. Its small size allows me to bring it with me to every appointment. In turn, I have been able to significantly increase the quality of care that I provide by diagnosing problems quickly and easily with certainty. ”

Clarius C7 Vet HD

- *Wireless Freedom*
- *High-Definition Imaging*
- *Easy App for iOS & Android*
- *Dedicated Small Animal Presets*
- *Affordable Scanner*
- *No Subscription Fees*
- *Free Clarius Cloud Storage*
- *Free Clarius Live Telemedicine*
- *Unlimited Users*



Questions?



Dr. Camilla Edwards



Dr. Oron Frenkel



Thank you!