

**WEBINAR**

# Visualize Your Needle: Ultrasound for Painless and Accurate MSK Injections

October 2022



# Your Host



## **Dr. Oron Frenkel, MD, MS**

Emergency Physician & POCUS Educator

Chairman, Clarius Medical Advisory Board

# Ultrasound-Guided Versus Blind Subacromial-Subdeltoid Bursa Injection

Review > Semin Arthritis Rheum. 2015 Dec;45(3):374-8.  
doi: 10.1016/j.semarthrit.2015.05.011. Epub 2015 May 21.

## Ultrasound-guided versus blind subacromial-subdeltoid bursa injection in adults with shoulder pain: A systematic review and meta-analysis

Tao Wu<sup>1</sup>, Hai Xin Song<sup>2</sup>, Yan Dong<sup>2</sup>, Jian Hua Li<sup>2</sup>

Affiliations + expand

PMID: 26590864 DOI: 10.1016/j.semarthrit.2015.05.011

### Abstract

**Objective:** This systematic review and meta-analysis aimed to assess the effectiveness of ultrasound-guided (USG) versus blind (landmark-guided, LMG) corticosteroid subacromial-subdeltoid bursa injection in adults with shoulder pain.

**Methods:** The searches were performed on PubMed, Ovid MEDLINE, Ovid EMBASE, Ovid CochraneCENTRAL, Web of Science, Google Scholar, and Scopus from database inception through March 27, 2015. Studies were included trials comparing USG versus LSG injections for the treatment of adults with subacromial-subdeltoid bursitis. Two reviewers independently performed data extraction and appraisal of the studies. The outcome measures collected were the decreased VAS and SDQ scores, the increased shoulder function scores and shoulder abduction motion range, and the effective rate at 6 weeks after injection.

**Results:** Seven papers including 445 patients were reviewed; 224 received LMG injections and 221 received USG injections. There was a significant difference in favor of USG for

*“Ultrasound-guided corticosteroid injections potentially offer a significantly greater clinical improvement over blind SASD bursitis injections in adults with shoulder pain.”*

Korbe S, Udoji EN, Ness TJ, Udoji MA. Ultrasound-guided interventional procedures for chronic pain management. Pain Manag. 2015;5(6):465-82. doi: 10.2217/pmt.15.46. Epub 2015 Sep 24. PMID: 26402316; PMCID: PMC4976830. Source: <https://pubmed.ncbi.nlm.nih.gov/26402316/>

# Ultrasound-Guided musculoskeletal interventional procedures around the shoulder

*“most **Ultrasound-guided** musculoskeletal interventional procedures around the shoulder produce **better results** in terms of accuracy and clinical efficacy than those performed in a blinded fashion.”*

Tortora S, Messina C, Gitto S, Chianca V, Serpi F, Gambino A, Pedone L, Carrafiello G, Sconfienza LM, Albano D. Ultrasound-guided musculoskeletal interventional procedures around the shoulder. J Ultrason. 2021 Jun 7;21(85):e162-e168. doi: 10.15557/JoU.2021.0026. Epub 2021 Jun 18. PMID: 34258042; PMCID: PMC8264815.



# Palpation Versus Ultrasound-Guided Acromioclavicular Joint Intra-articular Corticosteroid Injections: A Retrospective Comparative Clinical Study

*“US-guided AC joint IA injection for the treatment of symptomatic AC joint OA resulted in better pain and functional status improvement than palpation-guided IA injection.”*



Park KD, Kim TK, Lee J, Lee WY, Ahn JK, Park Y. Palpation Versus Ultrasound-Guided Acromioclavicular Joint Intra-articular Corticosteroid Injections: A Retrospective Comparative Clinical Study. Pain Physician. 2015 Jul-Aug;18(4):333-41. Erratum in: Pain Physician. 2015 Sep-Oct;18(5):517. PMID: 26218936.



# Ultrasound-Guided Suprapectoral Tenodesis of the Long Head of the Biceps Brachii

*“**Ultrasound** provides surgeons with a safe **and non-invasive tool** to **visualize the biceps tendon** as it exits the bicipital groove, negating the need for unroofing and other pitfalls associated with traditional techniques.”*


Andersen WJ, Barcelos M, de Paiva Raffaelli M, Hirahara AM. Ultrasound-Guided Suprapectoral Tenodesis of the Long Head of the Biceps Brachii. *Arthrosc Tech.* 2020 Dec 21;9(12):e2071-e2076. doi: 10.1016/j.eats.2020.08.039. PMID: 33381421; PMCID: PMC7768307.





# Poll

*Do you use ultrasound to guide procedures? If so, how often?*



# Your Expert Guest Speaker



**Alan M Hirahara, MD, FRCSC**

*Orthopaedic Surgery | Sports Medicine*



# EXPERT ULTRASOUND GUIDANCE FOR ACCURATE MSK INJECTIONS

## PART 1: THE SHOULDER

Alan M Hirahara, MD, FRCSC



SPORTS MEDICINE / SHOULDER & KNEE  
ARTHROSCOPY / ULTRASOUND / ORTHOBIOLOGICS



# Disclosures

- Consultant / Royalties
- Consultant
- Speaker / Stock options
- Committee member
  - CSMAS
  - Health & Safety
  - Industry
  - Development
- Arthrex Inc.
- LifeNet Health, Inc.
- Clarius Mobile Health
  - National Collegiate Athletic Association
  - Big Sky Athletic Conference
  - ASES Foundation
  - AANA



# Objectives

- How to accurately inject into the biceps and subscapularis tendons or sheaths
- Visualization of supraspinatus tendon abnormalities and PRP injection guidance
- The best approach for needle visualization during a subacromial space injection
- Accurate and effective AC joint and intra-articular shoulder injections
- Tips to pinpoint the optimal location for successful suprapectoral biceps tenodesis






# SHOULDER INJECTIONS



# Improved Accuracy


- Despite overwhelming evidence that accuracy in injections is quite low and ultrasound significantly improves this accuracy, many **DO NOT BELIEVE**
- Many physicians will feel that they are “excellent” at injections and “never” miss





# Shoulder Injections

	USG	Landmark guided
ACJ	93.6%	68.2%
Bicipital Groove	86.7%	26.7%
GHJ	92.5%	72.5%
Subacromial	65.0%	70.0% no stat significance




“Ultrasound-guided shoulder girdle injections are more accurate and more effective than landmark-guided injections: a systematic review and meta-analysis.”

Aly AR et al, Br J Sports Med. 2015




# My Experience

- “My patients have significantly less pain when done under ultrasound guidance”
- My Explanation:
  - ***“I was missing!”***



# Anisotropy



- **Sonographic artifact** associated with linear structures
  - Occurs when a linear structure is **not 90°** to ultrasound beam
  - Muscle, Tendon, Ligament, Cortex and Nerve
  - Tilting the transducer will cause the structure to appear artificially absent
  - **Shoulder Imaging PEARL**: If the **bone is bright**, the tendon typically is also







# Injections

In Plane



Out of Plane



# HOW DO I DO INJECTIONS?





# What Do You Need?

- Sterile gloves / gown / mask / hat
- Sterile drapes
- Sterile transducer sheath or cover
- Sterile gel
- Sterile operating room



# Sterile Technique




# BICEPS TENDON / GROOVE INJECTIONS






# Biceps Sheath Injection


Out of Plane



In Plane



# Biceps Tendon Injection







# SUPRASPINATUS / SUBACROMIAL INJECTIONS



# Subacromial Injection



# Supraspinatus / Subacromial Injections







# **GLENOHUMERAL JOINT INJECTIONS**



# Shoulder Joint Injection





# Glenohumeral Joint Injection



# AC JOINT INJECTIONS





# AC joint injection long axis needle





# AC Separation – Injection





# AC Joint Injection



# SURGICAL APPLICATIONS






# Ultrasound Assisted Surgery


- Take advantage of the tools:
  - Can identify anatomy / pathology
  - Can put a needle on it
  - Can use the centerline function
  - Can use the crosshairs

**–IN REAL TIME!**




# Ultrasound-Guided, Percutaneous Placement



In plane needle



Out of plane needle



Crosshair markings





# Publications

- A Guide to Ultrasound of the Shoulder, Part 3: Interventional & Procedural Uses
  - Hirahara A, Panero A
  - Am J Orthop, Nov-Dec 2016



## Orthopedic Technologies & Techniques

### A Guide to Ultrasound of the Shoulder, Part 3: Interventional and Procedural Uses

Alan M. Hirahara, MD, FRCS(C), and Alberto J. Panero, DO

#### Abstract

Ultrasound is an extremely useful diagnostic tool for physicians, but recent advances have found that ultrasound's greatest utility is in interventional and procedural uses. Numerous studies have demonstrated a significant improvement in outcome and patient satisfaction when using ultrasound guidance for injections. Newer techniques are emerging to use ultrasound as an aid to surgery and interventional procedures. This allows the physician to use smaller incisions and less invasive methods, which are also easier to use for the practitioner and more cost-effective.

for injections.<sup>2,22</sup> Within the limitation of using a needle, second-generation procedures—hydrodissection of peripherally entrapped nerves, capsular distention, mechanical disruption of neovascularization, and needle fenestration or barbotage in chronic tendinopathy—try to simulate surgical objectives while minimizing tissue burden and other complications of surgery.<sup>4</sup> More advanced procedures include needle fenestration/release of the carpal ligament in carpal tunnel syndrome and A1 pulley nodule release in the setting of trigger finger.<sup>5</sup> Innovative third-generation procedures involve the use of surgical tools such as hook blades under ultrasound guidance to perform surgical procedures. Surgeons are now improving already established percutaneous, arthroscopic, and open surgical procedures with ultrasound assistance.<sup>3</sup> Aside from better guidance, reducing cost and improving surgeon comfort may be additional benefits of ultrasound-assisted surgery.

Ultrasound has classically been marketed and used as a diagnostic tool. Radiologists, emergency physicians, and sports physicians used ultrasound units to rapidly and appropriately diagnose numerous injuries and disorders, in a timely and cost-effective manner. Part 1 and Part 2 of this series showed how to use ultrasound in the shoulder for diagnosis and how to code and get reimbursed for its use. Ultrasound can also be used to help guide procedures and interventions performed to treat patients. Currently, more physicians are beginning to recognize the utility of this modality as an aid to interventional procedures.

First-generation procedures use ultrasound to improve accuracy of joint, bursal, tendon, and muscular injections.<sup>2</sup> Recent studies have shown a significant improvement in accuracy, outcomes, and patient satisfaction using ultrasound guidance

#### Image-Guided Treatment Options

Prior to image guidance, palpation of surface anatomy helped physicians determine the anatomic placement of injections, incisions, or portals. Joints and bursas that do not have any inflammation or fluid can sometimes be difficult to identify or locate by palpation alone. Palpation-guided joint injections often miss their target and cause significant pain when the therapeutic agent is injected into a muscle, tendon, ligament, fat, or other tissue. Ultrasound-guided injections have proven to be more accurate and have better patient satisfaction when compared to blind injections.<sup>2,15</sup>

X-ray fluoroscopy has been the primary option for surgeons to assist in surgery. This is a natural modality for orthopedic surgeons; their primary use is for bone to help with fracture reduction and fixation as the bone, instrumentation, and

**Authors' Disclosure Statement:** Dr. Hirahara reports that he receives support from ArthroX as a consultant, royalties, and research support. Dr. Panero reports no actual or potential conflict of interest in relation to this article.

www.amjorthopedics.com

November/December 2016 The American Journal of Orthopedics® 1

# Publications

- Hirahara, Andersen. Ultrasound-guided percutaneous reconstruction of the *anterolateral ligament*: Surgical technique & case report. Am J Orthop. 2016 Nov-Dec;45(7):418-60.
- Hirahara, Andersen. Ultrasound-guided percutaneous repair of the *medial patellofemoral ligament*: Surgical technique & outcomes. Am J Orthop. May-June 2017;46(3):152-157.
- Hirahara A, Mackay G, Andersen W. “Ultrasound-guided suture tape augmentation and stabilization of the *medial collateral ligament*”, Arthrosc Tech, 2018 Feb 5; 7(3): e205-10. doi: 10.1016/j.eats.2017.08.069. PMID: 29881691
- Andersen J, Barcelos M, Raffaelli M, Hirahara A. “Ultrasound-guided suprapectoral tenodesis of the long head of the *biceps brachii*”, Arthrosc Tech, 2020 Dec 21; 9(12): e2071-6. doi: 10.1016/j.eats.2020.08.039. PMID: 33381421



# Ultrasound Assisted Shoulder Surgery

- Biceps tenodesis
  - Find the groove, place the portals
- Calcific tendonitis
- Loose bodies






# What's in My Bag?

- Linear, high-definition scanner
- Frequency: 5 – 15 MHz
- Depth Max: 7 cm
- Connection: Wireless (no cords)





# Suprapec Arthroscopic Biceps Tenodesis



## Technical Note

### Ultrasound-Guided Suprapectoral Tenodesis of the Long Head of the Biceps Brachii

Wyatt J. Andersen, M.S.H.S., A.T.C., Matheus Barcelos, M.D.,  
Maurício de Paiva Raffaelli, M.D., and Alan M. Hirahara, M.D., F.R.C.S.C.

**Abstract:** When the long head of the biceps tendon is diseased, tenodesis is an appropriate treatment strategy. The specific technique used is dependent on visualization, fixation method and hardware, and tenodesis location. For suprapectoral tenodesis techniques, those that fix the tendon within or below the bicipital groove can be challenging, owing to the transverse humeral ligament covering the groove. To accurately identify the biceps tendon in this area, the ligament often requires resection. Ultrasound provides surgeons with a safe and noninvasive tool to visualize the biceps tendon as it exits the bicipital groove, negating the need for unroofing and other pitfalls associated with traditional techniques. This technical note describes an ultrasound-guided suprapectoral biceps tenodesis procedure.

Numerous biceps tenodesis procedures have been described in the literature, and these procedures can vary based on the method of visualization, tenodesis location with respect to the bicipital groove and pectoralis major, and hardware.<sup>1-5</sup> Techniques that fix the long head of the biceps tendon (LHBT) at any position above the pectoralis major are typically performed arthroscopically, whereas a subpectoral tenodesis is an open or mini-open procedure.<sup>1,5</sup> For arthroscopic procedures performed in a suprapectoral manner, exposing the location in or below the groove can be challenging. The tissue layer overlying the bicipital groove (transverse humeral ligament) makes exposure of the tendon within the groove the greatest challenge using arthroscopy. For a suprapectoral tenodesis, the surgeon must visualize the biceps

tendon arthroscopically in the joint and tag the transverse humeral ligament at the top of the groove. In the subacromial space, this tag suture is used as the reference point to start unroofing the biceps, taking down the transverse humeral ligament and exposing the biceps within the groove.

Intraoperative use of ultrasound avoids these pitfalls. Ultrasound allows surgeons to accurately identify the LHBT as it exits the bicipital groove instead of going through the joint and unroofing the tissue layer above the groove. Ultrasound is an inexpensive, non-intrafascial, and noninvasive modality that can be easily used to visualize soft-tissue structures in real time. In this technical note, we describe ultrasound-guided suprapectoral tenodesis of the LHBT (Video 1, Table 1).

*From Private Practice, Sacramento, California, U.S.A. (W.J.A., A.M.H.) and Institute of the Shoulder and the Elbow, Nílkro Assunção de Emílio em Ortopedia e Neurocirurgia (P.A.R.), São Paulo, Brazil (M.B., M.A.P., F.J.). The authors report the following potential conflicts of interest or sources of funding: A.M.H. reports royalties from Arthrosc; is a consultant for Arthrosc and Ligoon Health; receives research support from Arthrosc; is a medical advisor for Clarus Medical Health; and owns stock options in Clarus Medical Health, outside the submitted work. In addition, A.M.H. has a patent with Arthrosc (United States) pending for TENOSURABLE CONSTRUCTS WITH MULTILIMB LOCKING MECHANISM THROUGH SINGLE SPINCE AND METHODS OF TISSUE REPAIR; has a patent with Arthrosc (European Patent Commission) pending for SUTURE ANCHOR AND METHODS OF KNITLESS TISSUE RECONSTRUCTION TECHNIQUES; has a patent with Arthrosc pending for TENOSURABLE CONSTRUCTS WITH MULTILIMB LOCKING MECHANISM THROUGH SINGLE SPINCE; has a patent with Arthrosc (European Patent Commission) pending for TENOSURABLE CONSTRUCTS WITH MULTILIMB LOCKING MECHANISM THROUGH SINGLE*

*SPINCE; has patents with Arthrosc (United States) issued and issued for SUTURE ANCHOR AND METHODS OF KNITLESS TISSUE FIXATION; has patents with Arthrosc (United States) issued and issued for MEASURING TOOL, USING SUTURE AND SUTURE ANCHOR; and has patent with Arthrosc (United States) pending for JOINT KINEMATIC RECONSTRUCTION TECHNIQUES. Full disclosures and disclosures forms are available for this article online, as supplementary material. Received June 22, 2020; accepted August 29, 2020. Address correspondence to Alan M. Hirahara, M.D., F.R.C.S.C., Private Practice, 2801 K St, SW 300, Sacramento, CA 95816, U.S.A. E-mail: alanhirahara@arthrosc.com. © 2020 by the Arthroscopy Association of North America. Published by Elsevier. This is an open access article under the CC BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/>). 2212-6287/2010/\$ - see front matter © 2020 AAOS. <https://doi.org/10.1016/j.atsk.2020.08.019>*



# Ultrasound-Guided Biceps Tenodesis



# Summary of Ultrasound

- *Tool* like fluoroscopy or arthroscopy
- *Inexpensive* and widely available
- Allows for *visualization* of anatomic structures
- Helps *guide* injections or surgeries
- Transforms surgeries from open to *percutaneous*





THANK YOU!



SPORTS MEDICINE / SHOULDER & KNEE  
ARTHROSCOPY / ULTRASOUND / ORTHOBIOLOGICS



# “Live Demonstration



**Shelley Guenther**  
*Clinical Marketing Manager*



# *Poll*

*What additional information would you like?*





# Clarius I



Genese Castonguay

- *30% Smaller & More Affordable*
- *Wireless Freedom*
- *High-Definition Imaging*
- *Easy App for iOS & Android*
- *Needle Enhance*
- *Clarius Cloud Storage*
- *Clarius Live Telemedicine*
- *Unlimited Users*



# ***Poll: Pre-Register***

***FREE WEBINAR***

***Expert Ultrasound Guidance  
for Accurate MSK Injections,  
Part 2: The Knee***

*January 24th, 2022*

*2 PM Pacific | 5 PM Eastern*

[www.clarius.com/ultrasound-webinars](http://www.clarius.com/ultrasound-webinars)

# Questions



*Dr. Hirahara*



*Dr. Frenkel*





*Thank you!*