

WEBINAR

Ultrasound-Guided Techniques to Rapidly Dissolve Fillers, Part 1: Misplaced Facial Fillers

January 2023



Your Host



Shelley Guenther, CRGS, CRCS

Clinical Marketing Manager

Ultrasound to Improve the Safety of Hyaluronic Acid Filler Treatments

> J Cosmet Dermatol. 2018 Dec;17(6):1019-1024. doi: 10.1111/jocd.12726. Epub 2018 Aug 6.

Ultrasound to improve the safety of hyaluronic acid filler treatments

Leonie W Schelke¹, Tom S Decates¹, Peter J Velthuis¹

Affiliations + expand
PMID: 30084182 DOI: 10.1111/jocd.12726

Abstract

Background: Hyaluronic acid fillers are known for a reliable safety profile, but complications do occur, even serious vascular adverse events.

Objective: To improve the safety of hyaluronic acid filler treatments.

Methods: Ultrasound is used to image hyaluronic acid fillers.

Results: Before a filler treatment is performed with ultrasound, previous filler treatments can be brought in to sight and vascular mapping can be performed. In case of adverse events, the filler and the surrounding tissues are visible. Dislocation, abscesses, and vascular adverse events can be seen. Under ultrasound guidance, hyaluronidase can be injected directly into the filler deposit.

Conclusion: Ultrasound examination can be an important tool to improve the safety of hyaluronic acid filler treatments.

Keywords: complications; cosmetic dermatology; filler; hyaluronic acid; safety; ultrasound.

J Cosmet Dermatol. Published by Wiley Periodicals, Inc.

“Dislocation, abscesses, and vascular adverse events can be seen. Under ultrasound guidance, hyaluronidase can be injected directly into the filler deposit.”

Schelke LW, Decates TS, Velthuis PJ. Ultrasound to improve the safety of hyaluronic acid filler treatments. J Cosmet Dermatol. 2018 Dec;17(6):1019-1024. doi: 10.1111/jocd.12726. Epub 2018 Aug 6. PMID: 30084182.. Source: <https://pubmed.ncbi.nlm.nih.gov/30084182/>

Ultrasound patterns of different dermal filler materials used in aesthetics

> J Cosmet Dermatol. 2021 May;20(5):1541-1548. doi: 10.1111/jocd.14032. Epub 2021 Mar 11.

Ultrasound patterns of different dermal filler materials used in aesthetics

Fernando Urdiales-Gálvez ¹, Francisco M De Cabo-Francés ², Isabel Bové ²

Affiliations + expand

PMID: 33641224 PMCID: PMC8252486 DOI: 10.1111/jocd.14032

[Free PMC article](#)

Abstract

Background Hyaluronic acid (HA) injection procedures has experienced an unprecedented increase. Aims To assess and determine, by using ultrasound examinations, the patterns corresponding to different dermal fillers. Patients/Methods Observational and retrospective bicenter study conducted on patients who underwent previous aesthetic treatments with dermal fillers. Ultrasound examinations were performed, at each study center, by one experienced observer. Results Sixty patients were included in the analysis. Among them, 48 patients showed a well-defined ultrasound pattern, while 12 exhibited a mixed one. According to ultrasound images, 4 different patterns were identified: [1] Heterogeneous, characterized by alternating hyperechoic and anechoic areas, which are visualized in the tissue in a heterogeneous way. This pattern is associated with healthy skin/subcutaneous cellular tissue and with fully integrated HA fillers. [2] Fine grain snowfall, characterized by alternating hyperechoic imaging, with posterior echogenic shadows. It is typical of liquid injectable silicone. [3] Coarse grain snowfall, characterized by

*"The identification of these patterns will allow specialists to choose the **best therapeutic approach** in patients who underwent **previous aesthetic treatments**"*

Urdiales-Gálvez F, De Cabo-Francés FM, Bové I. Ultrasound patterns of different dermal filler materials used in aesthetics. J Cosmet Dermatol. 2021 May;20(5):1541-1548. doi: 10.1111/jocd.14032. Epub 2021 Mar 11. PMID: 33641224; PMCID: PMC8252486.

Filler Migration: A Number of Mechanisms to Consider

"It is important all physicians assessing nodules/masses/swelling in the facial area be aware that soft tissue fillers may migrate to a location away from their intended site of injection by several mechanisms and persist in the tissue even years later."

Jordan DR, Stoica B. Filler Migration: A Number of Mechanisms to Consider. Ophthalmic Plast Reconstr Surg. 2015 Jul-Aug;31(4):257-62. doi: 10.1097/IOP.0000000000000368. PMID: 25650796.

Filler Migration: A Number of Mechanisms to Consider

David R Jordan¹, Bazil Stoica

Affiliations + expand

PMID: 25650796 DOI: [10.1097/IOP.0000000000000368](https://doi.org/10.1097/IOP.0000000000000368)

Abstract

Purpose: To report 3 representative cases of soft tissue filler identified in locations other than their intended injected sites (possible migration) and review the literature on pathogenesis of filler migration.

Introduction: Soft tissue fillers are continuing to increase in popularity throughout North America and worldwide as a means of volume restoration and contour enhancement. With increasing recognition of their value in restoring a more youthful appearance and the ease of office injection, soft tissue fillers have become one of the most commonly performed nonsurgical cosmetic procedures. Soft tissue fillers are also foreign bodies in our system and therefore have the potential for a myriad of complications both immediately after the injection and potentially months or years later. Filler migration is one such complication and has a number of potential mechanisms.

Methods: The authors reviewed the medical records of 3 patients with filler located in areas other than their intended injected sites possibly as a result of migration. All patients were from the practice of 1 individual (DRJ). A MEDLINE search of the English-language literature on filler migration was conducted to investigate the various causes responsible for migration of filler.

Results: Clinical manifestations of the possible filler migration in the 3 cases included eyelid swelling in 2 patients and a noninflammatory mass adjacent to the area of filler injection in the third patient. Surgery was performed on 1 patient, and filler was visualized in the tissue and dissolved with hyaluronidase. Hyaluronidase was also used to dissolve the suspected filler in the third patient, and the third patient has elected to continue with nonsurgical management.

Conclusions: Filler migration is a potential complication of soft tissue filler injection. Physicians should be aware of the various mechanisms of filler migration and the potential for complications. Filler migration should be considered in the differential diagnosis of nodules, masses, or swelling in the facial area.

European Federation of Societies for Ultrasound in Medicine Position Statement on Dermatologic Ultrasound

*“a **position statement** from a scientific society with regards to this **application of US** is useful for physicians involved in dermatologic US, allowing for the foundations of present clinical practice and **evidence generation** to be developed.”*

Fernando Alfageme¹, Ximena Wortsman², Orlando Catalano³, Gaston Roustan¹, Maria Crisan⁴, Diana Crisan⁵, Diana E. Gaitini⁶, Eugenio Cerezo⁷, Radu Badea⁸



European Federation of Societies for Ultrasound in Medicine Position Statement on Dermatologic Ultrasound

“ POSITION STATEMENT 2 The minimum transducer frequency for dermatologic ultrasound should be 15 MHz. Higher transducer frequencies may provide further information that may be relevant (broad agreement 8/9, 88.9 %)..”

Fernando Alfageme¹, Ximena Wortsman², Orlando Catalano³, Gaston Roustan¹, Maria Crisan⁴, Diana Crisan⁵, Diana E. Gaitini⁶, Eugenio Cerezo⁷, Radu Badea⁸



European Federation of Societies for Ultrasound in Medicine Position Statement on Dermatologic Ultrasound

“US can provide relevant information that includes data on facial anatomical variants, the type, location and extent of common cosmetic fillers, the identification of implants, the complications of lipolytic procedures, and the possibility of percutaneous US guidance for the procedure .”


Fernando Alfageme¹, Ximena Wortsman², Orlando Catalano³, Gaston Roustan¹, Maria Crisan⁴, Diana Crisan⁵, Diana E. Gaitini⁶, Eugenio Cerezo⁷, Radu Badea⁸





Poll

What key benefits do you see ultrasound bringing to facial aesthetics?



Your Expert Guest Speaker



Dr. MJ Rowland-Warmann

*BDS, BSc, MSc, Aes. Med, MClin, Dent Orthod,
PGDip Endod, PGCert MedEd, Dip.MJDF*

Founder Smileworks Liverpool

Ultrasound guided techniques to dissolve fillers, Part 1: Misplaced facial fillers

Dr MJ Rowland-Warmann

BDS BSc MSc Aes.Med. MClinDent Orthod. PGCert MedEd. Dip.MJDF

Founder & Lead Clinician Smileworks HUB



dr_mj_smileworks



L20

How ultrasound can help us

- Filler treatment – before, during and after
- Monitoring filler over time
- Diagnosing filler types
- Learning anatomy
- Diagnosis and management of complications

How ultrasound can help us

- Filler treatment – before, during and after
- Monitoring filler over time
- Diagnosing filler types
- Learning anatomy
- **Diagnosis and management of complications**

Misplaced filler

- Complications
 - Immediate
 - Delayed
- At the time of treatment & filler spread
- Common
 - patient impact
 - Industry impact



Funt, D. and T. Pavicic (2013). "Dermal fillers in aesthetics: an overview of adverse events and treatment approaches." Clin Cosmet Investig Dermatol 6: 295-316.

Early events (occurring up to several days post-treatment)	Delayed events (occurring from weeks to years post-treatment)
Injection site reactions <ul style="list-style-type: none"> Erythema Edema Pain/tenderness Bruising Itching 	Infection (atypical; eg, mycobacterial) <ul style="list-style-type: none"> Erythema Edema Pain/tenderness Nodule/abscess Systemic responses to infection Biofilm
Infection <ul style="list-style-type: none"> Erythema Edema Pain/tenderness Acne papule formation Nodule/abscess 	Foreign body granuloma <ul style="list-style-type: none"> Varying from subclinical histologic changes to disfiguring nodules
Hypersensitivity <ul style="list-style-type: none"> Erythema Edema Pain/tenderness Nonfluctuant nodules 	Migration of implant material
Lumps, asymmetries, contour irregularities caused by technique and placement errors	Immune reactions <ul style="list-style-type: none"> Local and site of injection and generalized Persistent discoloration Persistent scarring
Skin discoloration <ul style="list-style-type: none"> Redness Whiteness Hyperpigmentation 	
Local tissue necrosis caused by vascular occlusion	Malar edema


Tear trough



Tear trough: layered anatomy






Tear trough: vascular anatomy



Anido, J., et al. (2021).
"Recommendations for the treatment of
tear trough deformity with cross-linked
hyaluronic acid filler." *J Cosmet
Dermatol* **20**(1): 6-17.



Techniques for tear trough



Tear trough complications

- Delayed
 - Swelling (42.3%)
 - Nodules/lumps (25%)
 - Migration (7.7%)

Delayed Complications following Dermal Filler for Tear Trough Augmentation: A Systematic Review

Lily Nguyen Trinh, MD^{1,2}  Kelly C. McGuigan, MS³ Amar Gupta, MD⁴ 

¹ School of Medicine, Tulane University School of Medicine, New Orleans, Louisiana


² Department of Otolaryngology—Head and Neck Surgery, Mass Eye and Ear, Boston, Massachusetts

³ School of Medicine, Thomas Jefferson University Sidney Kimmel Medical College, Philadelphia, Pennsylvania

⁴ Department of Otolaryngology, Private Practice—Head and Neck Surgery, Los Angeles, California

Address for correspondence Lily Nguyen Trinh, MD, Tulane University School of Medicine—School of Medicine, 1430 Tulane Avenue, New Orleans, LA 70112-2632 (e-mail: lilytrinh4@gmail.com).

Pre-operative imaging




Video

Guided hyaluronidase



Video

Immediately after





Video

Forehead



Forehead: vascular anatomy

Beleznay, K., et al. (2019).
"Update on Avoiding and
Treating Blindness From Fillers:
A Recent Review of the World
Literature." *Aesthet Surg J*
39(6): 662-674.



Forehead: layered anatomy



Pre-operative scanning



Video

Guided hyaluronidase



Video

1 hour after hyaluronidase



Video

Complication management

- Incorporate ultrasound into your daily routine
- Build your skills gradually
- Start with “everyday” complications
- Ask for help if you are not sure




FACIAL ULTRASOUND TRAINING

use ultrasound to plan treatments, reduce complications and inject safely



START HERE



20% off online course only. Valid for 30 days:

CLARIUSWEBINAR20

Thank you



@dr_mj_smileworks
@smileworkshub

Email: delegates@smileworks-hub.co.uk

Web: www.smileworks-hub.co.uk

≡ **Aesthetic HUB**
Liverpool, UK



“Live Demonstration



Shelley Guenther

Clinical Marketing Manager

First Handheld Ultrasound with Limitless Scanning

 clarius





Clarius L20 HD₃

30% Smaller & More Affordable
Wireless Freedom

High-Definition Imaging

Easy App for iOS & Android

Limitless Scanning with Power Fan
HD₃

Advanced Aesthetics Package

Clarius Cloud Storage

Clarius Live Telemedicine

Unlimited Users

Questions?



*Dr. MJ Rowland-
Warmann*



Shelley Guenther



Thank you!