



WEBINAR

Pragmatic MSK Ultrasound: 8-Step Protocol for Assessing the Distal Biceps Tendon

October 2023



Your Host



Shelley Guenther, CRGS, CRCS

Sonographer | Clinical Marketing Manager

Ultrasound assessment in distal biceps tendon injuries: Techniques, pearls and pitfalls

“

Ultrasound of the distal BT is challenging.... a combination of approaches and examination techniques are utilized to ensure complete visualization.”

Al-Ani Z, Lauder J. Ultrasound assessment in distal biceps tendon injuries: Techniques, pearls and pitfalls. Clin Imaging. 2021 Jul;75:46-54. doi: 10.1016/j.clinimag.2021.01.017. Epub 2021 Jan 20. PMID: 33493736.

Review > Clin Imaging. 2021 Jul;75:46-54. doi: 10.1016/j.clinimag.2021.01.017.

Epub 2021 Jan 20.

Ultrasound assessment in distal biceps tendon injuries: Techniques, pearls and pitfalls

Zeid Al-Ani¹, Joshua Lauder²

Affiliations + expand

PMID: 33493736 DOI: [10.1016/j.clinimag.2021.01.017](https://doi.org/10.1016/j.clinimag.2021.01.017)

Abstract

Ultrasound assessment of the distal biceps tendon is challenging. The tendon has two components which are continuations of the long and short heads of the muscle, and these undergo 90° of rotation along their course. The tendon has a deep insertion to the radial tuberosity. Therefore, a combination of approaches and examination techniques are utilized to ensure complete evaluation. The various ultrasound approaches used to assess the distal biceps tendon, with their advantages and limitations, will be described. Selected examples of distal biceps tendon injuries, with magnetic resonance imaging (MRI) correlation in challenging cases, will be demonstrated.

Keywords: Biceps Brachii; Distal biceps tendon; Muscle; Skeletal; Tendon injuries; Ultrasonography.

Copyright © 2021 Elsevier Inc. All rights reserved.

[PubMed Disclaimer](#)

Similar articles

[Partial rupture of the distal biceps brachii tendon: a magnetic resonance imaging analysis.](#)

Nicolay RW, Lawton CD, Selley RS, Johnson DJ, Vassa RR, Prescott AE, Omar IM, Marra G.

J Shoulder Elbow Surg. 2020 Sep;29(9):1859-1868. doi: 10.1016/j.jse.2020.04.021. Epub 2020 Jun 9.

PMID: 32815807

[Terminal bifurcation of the biceps brachii muscle and tendon: anatomic considerations and clinical implications.](#)

Dirim B, Brouha SS, Pretterklieber ML, Wolff KS, Frank A, Pathria MN, Chung CB.

AJR Am J Roentgenol. 2008 Dec;191(6):W248-55. doi: 10.2214/AJR.08.1048.

PMID: 19020211

[Distal Biceps Brachii Tendon Insertion: A Simple Method of Ultrasound Evaluation.](#)

Draghi F, Bortolotto C, Ferrozzi G.

J Ultrasound Med. 2021 Apr;40(4):811-813. doi: 10.1002/jum.15459. Epub 2020 Aug 31.

Distal Biceps Brachii Tendon Insertion: A Simple Method of Ultrasound Evaluation

“

The **crab position**... allows for **optimal visualization** of the **distal biceps tendon** in the axial plane ”

Draghi F, Bortolotto C, Ferrozzi G. Distal Biceps Brachii Tendon Insertion: A Simple Method of Ultrasound Evaluation. J Ultrasound Med. 2021 Apr;40(4):811-813. doi: 10.1002/jum.15459. Epub 2020 Aug 31. PMID: 32865288.

> J Ultrasound Med. 2021 Apr;40(4):811-813. doi: 10.1002/jum.15459. Epub 2020 Aug 31.

Distal Biceps Brachii Tendon Insertion: A Simple Method of Ultrasound Evaluation

Ferdinando Draghi ^{1 2}, Chandra Bortolotto ², Guia Ferrozzi ^{3 4}

Affiliations + expand

PMID: 32865288 DOI: 10.1002/jum.15459

Abstract

Different ultrasound approaches have been proposed to evaluate the distal biceps brachii tendon, such as lateral, posterior, medial, and anterior, all allowing good visualization of the distal biceps tendon and its pathologic appearances. Here we describe a new method of ultrasound evaluation of the distal biceps tendon insertion with the patient in the so-called crab position. The crab position is used to evaluate the posterior compartment of the elbow but also the lateral and medial compartments. This position allows for coverage of 75% the elbow in a single position, with less discomfort for the patient and a fair amount of saved time. In this position, after the evaluation of the common extensor tendon in the long axis, turning the transducer 90° and moving it distally allows for optimal visualization of the distal biceps brachii tendon in the axial plane. It is a simple, quick approach that allows for excellent visualization of the insertion of the distal biceps brachii tendon while avoiding an anisotropic effect.

Keywords: distal biceps brachii tendon; tendon; ultrasound.

© 2020 American Institute of Ultrasound in Medicine.

[PubMed Disclaimer](#)

Similar articles

[Disorders of the distal biceps brachii tendon.](#)

Chew ML, Giuffrè BM.

Radiographics. 2005 Sep-Oct;25(5):1227-37. doi: 10.1148/rg.255045160.

PMID: 16160108 [Review](#).

[Ultrasound assessment in distal biceps tendon injuries: Techniques, pearls and pitfalls.](#)

Al-Ani Z, Lauder J.

Clin Imaging. 2021 Jul;75:46-54. doi: 10.1016/j.clinimag.2021.01.017. Epub 2021 Jan 20.

PMID: 33493736 [Review](#).

[Terminal bifurcation of the biceps brachii muscle and tendon: anatomic considerations and clinical implications.](#)

Dirim B, Brouha SS, Pretterklieber ML, Wolff KS, Frank A, Pathria MN, Chung CB.

Ultrasound classification of traumatic distal biceps brachii tendon injuries

“

Ultrasound can be used as a **first-line imaging modality** to evaluate distal biceps tendon injuries.”

de la Fuente J, Blasi M, Martínez S, Barceló P, Cachán C, Miguel M, Pedret C. Ultrasound classification of traumatic distal biceps brachii tendon injuries. *Skeletal Radiol.* 2018 Apr;47(4):519-532. doi: 10.1007/s00256-017-2816-1. Epub 2017 Nov 24. PMID: 29177701; PMCID: PMC5814508.

> [Skeletal Radiol.](#) 2018 Apr;47(4):519-532. doi: 10.1007/s00256-017-2816-1. Epub 2017 Nov 24.

Ultrasound classification of traumatic distal biceps brachii tendon injuries

Javier de la Fuente ¹, Marc Blasi ^{2 3}, Sílvia Martínez ⁴, Pablo Barceló ^{5 6}, Carlos Cachán ¹, Maribel Miguel ⁷, Carles Pedret ^{8 9}

Affiliations + expand

PMID: 29177701 PMCID: [PMC5814508](#) DOI: [10.1007/s00256-017-2816-1](#)

[Free PMC article](#)

Abstract

Objective: The present work is aimed at analysing ultrasound findings in patients with distal biceps brachii tendon (DBBT) injuries to assess the sensitivity of ultrasound in detecting the different forms of injury, and to compare ultrasound results with magnetic resonance imaging (MRI) and surgical results.

Materials and methods: A total of 120 patients with traumatic DBBT injuries examined between 2011 and 2015 were analysed. We compared ultrasound results with MRI results when surgery was not indicated and with MRI and surgical results when surgery was indicated.

Results: For major DBBT injuries (complete tears and high-grade partial tears), the concordance study between exploration methods and surgical results found that ultrasound presented a slight statistically significant advantage over MRI (ultrasound: $\kappa = 0.95$ -very good-95% CI 0.88 to 1.01, MRI: $\kappa = 0.63$ -good-95% CI 0.42 to 0.84, kappa difference $p < 0.01$). Minor injuries, in which most tendon fibres remain intact (tendinopathies, elongations and low-grade partial tears), are the most difficult to interpret, as ultrasound and MRI reports disagreed in 12 out of 39 cases and no surgical confirmation could be obtained.

Conclusions: Based on present results and previous MRI classifications, we establish a traumatic DBBT injury ultrasound classification. The sensitivity and ultrasound-surgery correlation results in the diagnosis of major DBBT injuries obtained in the present study support the recommendation that ultrasound can be used as a first-line imaging modality to evaluate DBBT injuries.

Keywords: Distal biceps tendon; External bicipital aponeurosis; Imaging technique; Injury classification; Tendon retraction; Ultrasonography.

[PubMed Disclaimer](#)


Conflict of interest statement

All the authors declare that they have neither conflicts of interest nor any disclosure to declare.



Interactive Poll

When you come across MSK injuries or abnormalities, Which of these challenges have you encountered in treating MSK concerns?



Your Expert Speaker



Marc Schmitz, MSc

Founder & CEO of Sonoskills

A background image showing a person's arm being examined with an ultrasound probe. The probe is held against the elbow area, and the arm is resting on a white surface. The text is overlaid on this image.

The 8-step distal biceps tendon ultrasound protocol

Marc Schmitz, MSc.

Disclosure

I have NO financial relationship or conflicts of interest with the presented material in this activity



Founder / Trainer




marc@sonoskills.com

[Linkedin.com/in/marcsonoskills](https://www.linkedin.com/in/marcsonoskills)


+7.000 ultrasound pathology cases



Regular contributor



free access!!











- Flat tendon
- No tendon sheath (and thus no tenosynovitis)
- One or two (or more) separate tendons (short/long head):
 - Myotendinous junction: both positioned side by side; long head lateral
 - Distal: changing orientation, demonstrating 90° of tendon external rotation; short head more superficial.
- Insertion:
 - Long head: proximally, with an oval shaped footprint to most of radial tuberosity
 - Short head: fanshaped and more distally into the radial tuberosity

- Size whole footprint:
 - Length: ± 21 mm
 - Width: ± 7 mm
- CSA:
 - Long head: ± 7.2 mm²
 - Short head: ± 5.2 mm²
- Lacertus fibroses (bicipital aponeurosis):
 - Stabilizes distal tendon
 - Load distribution
- Bicipital bursa:
 - Surrounds the attachment of both tendons
 - Acts to reduce friction between the tendons and the radial tuberosity

1. Clinical examination
2. Ventral DBT + BA approach: SAX
3. Ventral DBT approach: LAX
4. Dorsal DBT approach: LAX + dynamical test
5. Dorsal DBT approach: SAX
6. Medial DBT approach: LAX
7. Lateral DBT approach: LAX + dynamical test
8. Pathology classification ^{*SonoSkills pathology checklist}

Step 1

Clinical Examination

Reversed "Popeye" sign




Distal biceps provocation test










Diagnosis of complete DBT tears is **mainly clinical**, whereas in partial tears medical imaging is a valuable addition to the clinical diagnosis.

Step 2

Ventral approach of DBT and BA (SAX)

Step 3

Ventral approach of DBT (LAX)




Step 4

Dorsal approach of DBT + dynamical test (LAX)

Step 5

Dorsal approach of DBT (SAX)



Step 6

Medial approach of DBT (LAX)



Step 7

Lateral approach of DBT + dynamical test (LAX)



MUSCULOSKELETAL ULTRASOUND

Hands-on and online ultrasound education for Sonoholics

COURSES



Distal biceps tendon

(Sono)pathology

- Distal biceps tendon < proximal biceps tendon
- Rupture incidence: 1.2 per 100.000 per year
- Mostly: men between ages of 40 – 60 years
- Risk factors:
 - Weightlifting
 - Anabolic steroids
 - Smoking
 - Poor physical condition
 - Increasing activity levels men in 5th and 6th decades
 - Systemic processes:
 - Several rheumatology diseases
 - Chronic uremia
 - Hemodialysis
 - Hyperparathyroidism

- Complete tear as a result of single traumatic event:
 - Or > sudden eccentric load applied to a flexed elbow
 - Or > forceful hyperextension against resistance
- Sensation:
 - Pop / sharp pain in antecubital fossa
 - Weakness elbow flexion / supination
- Tears frequently result in functional weakness and disability
- Bruising
- 2 weeks: within this time preferred diagnosis and surgical repair
- Palpable defect, or tendon retraction, only if:
 - Distal biceps tendon is avulsed from radial tuberosity
 - AND torn lacertus fibrosis
- Site: mostly within 1 – 2 cm of its insertion
- ± 2 cm: location of hypovascular zone

- Partial tear
 - Uncommon: partial tear < complete tears
 - Underlying process maybe bone impingement or tendinosis
 - Single-bundle ruptures are presumably underreported
 - Difficult to diagnose:
 - Antecubital fossa pain
 - May or may not have a history of a traumatic event
 - A painful but intact tendon, verified by a hook test, is indicative of a partial tear
 - Some patients with bicipitoradial bursitis can present with an enlarging mass

- Pitfall 1: intact lacertus fibrosis makes clinical diagnosis more difficult.
 - Retraction <8 cm typically correlates with an intact lacertus
 - Retraction >8 cm indicates a torn lacertus
- Pitfall 2: if the distal tendon is bifurcated an isolated rupture of a single tendon, with the other remaining intact, is possible

- Mechanical enthesopathy
- Inflammatory enthesitis
- Calcifications in the tendon or the enthesis
- Neovascularization, e.g. by means of B-flow examination
- Cubital bursitis (bicipitoradial bursitis)
- Interosseous bursitis
- Tendinosis (usually 1–2 cm proximal)
- Partial ruptures
- Complete ruptures
- Avulsions with osseous tearing
- Postoperative assessment with evaluation of tendon integrity, anchor position, suture position
- Determination of tendon elasticity using elastography
- Mechanical impingement during pronation-supination
- Secondary irritation caused by enthesophytes, osteophytes or by articular pathologies (synovitis, activated cubital arthrosis, synovial chondromatosis, etc.)

Clinical indications for musculoskeletal ultrasound updated in 2017 by European Society of Musculoskeletal Radiology (ESSR) consensus

Pathology

Consensus

Distal biceps insertion

2

Bicipitoradial bursitis

2

Step 8

Pathology classification *SonoSkills pathology checklist

1. STRUCTURE OF PATHOLOGY



Shape – *did the morphology change?*



Echogenicity – *did the echogenicity of the anatomical structure change?*



Continuity – *are there any tears or fiber discontinuities?*



Doppler – *is there a presence of neovascularization?*



Functional – *how does the anatomy move?*

2. INTERPRET PATHOLOGY



LAX/SAX – *scan in two planes*



Left/right – *compare the affected with the non-affected side*



Measure – *quantify morphologic changes*




Relevance – *how clinical relevant are the changes?*

MRI classification system	Ultrasound classification system	Ultrasound description	Number of cases	Relative (%)
No tear	Type 1	Thickened, hypoechogenic, loss of the fibrillary pattern, tendon continuity to the attachment retained, no evidence of torn fibres on the static and dynamic examinations	20	16.66
	Type 1a (Fig. 7)	It particularly affects one of the two DBBT components (short head or long head)	4	3.33
	Type 1b (Fig. 8)	It affects the full thickness of both DBBT components (short head and long head)	16	13.33
Partial tear	Type 2	DBBT thinning, irregular contour, anechoic appearance, partial tendon discontinuity, peritendinous effusion	35	29.16
Low grade (≤ 50 tear)	Type 2a (Fig. 9)	Tendon discontinuity at the attachment site, thinning $\leq 50\%$ of the total thickness	18	15.00
High grade (>50 tear)	Type 2b (Fig. 5)	Tendon discontinuity at the attachment site, $>50\%$ of the total thickness	14	11.66
	Type 2c (Fig. 6)	Tendon discontinuity at the attachment site that affects the full thickness of a single DBBT component (short head or long head), which becomes retracted. A refraction artefact at the level of the tendinous stump is observed	3	2.50
Complete tear	Type 3	Proximal hypertrophic tendon stump with a refraction artefact, snake-like pattern on the long-axis view, absence of tendon fibres at the attachment site, peritendinous effusion	65	54.16
Non-retracted (≤ 8 cm)	Type 3a (Fig. 3)	EBA intact or hypertrophied, with minor tendon retraction	10	8.33
Retracted (>8 cm)	Type 3b (Fig. 4)	EBA tear with marked DBBT retraction, no evidence of EBA continuity	55	45.83


MRI classification system	Ultrasound classification system	Ultrasound description	Number of cases	Relative (%)
No tear	Type 1	Thickened, hypoechogenic, loss of the fibrillary pattern, tendon continuity to the attachment retained, no evidence of torn fibres on the static and dynamic examinations	20	16.66
	Type 1a (Fig. 7)	It particularly affects one of the two DBBT components (short head or long head)	4	3.33
	Type 1b (Fig. 8)	It affects the full thickness of both DBBT components (short head and long head)	16	13.33
Partial tear	Type 2	DBBT thinning, irregular contour, anechoic appearance, partial tendon discontinuity, peritendinous effusion	35	29.16
	Type 2a (Fig. 9)	Tendon discontinuity at the attachment site, thinning $\leq 50\%$ of the total thickness	18	15.00
	Type 2b (Fig. 5)	Tendon discontinuity at the attachment site, $>50\%$ of the total thickness	14	11.66
Complete tear	Type 2c (Fig. 6)	Tendon discontinuity at the attachment site that affects the full thickness of a single DBBT component (short head or long head), which becomes retracted. A refraction artefact at the level of the tendinous stump is observed	3	2.50
	Type 3	Proximal hypertrophic tendon stump with a refraction artefact, snake-like pattern on the long-axis view, absence of tendon fibres at the attachment site, peritendinous effusion	65	54.16
Non-retracted (≤ 8 cm)	Type 3a (Fig. 3)	EBA intact or hypertrophied, with minor tendon retraction	10	8.33
Retracted (>8 cm)	Type 3b (Fig. 4)	EBA tear with marked DBBT retraction, no evidence of EBA continuity	55	45.83

Distal biceps: Type 1a pathology

LAX



SAX




A





B

MRI classification system	Ultrasound classification system	Ultrasound description	Number of cases	Relative (%)
No tear	Type 1	Thickened, hypoechogenic, loss of the fibrillary pattern, tendon continuity to the attachment retained, no evidence of torn fibres on the static and dynamic examinations	20	16.66
	Type 1a (Fig. 7)	It particularly affects one of the two DBBT components (short head or long head)	4	3.33
	Type 1b (Fig. 8)	It affects the full thickness of both DBBT components (short head and long head)	16	13.33
Partial tear	Type 2	DBBT thinning, irregular contour, anechoic appearance, partial tendon discontinuity, peritendinous effusion	55	29.16
Low grade (≤ 50 tear)	Type 2a (Fig. 9)	Tendon discontinuity at the attachment site, thinning $\leq 50\%$ of the total thickness	18	15.00
High grade (>50 tear)	Type 2b (Fig. 5)	Tendon discontinuity at the attachment site, $>50\%$ of the total thickness	14	11.66
	Type 2c (Fig. 6)	Tendon discontinuity at the attachment site that affects the full thickness of a single DBBT component (short head or long head), which becomes retracted. A refraction artefact at the level of the tendinous stump is observed	3	2.50
Complete tear	Type 3	Proximal hypertrophic tendon stump with a refraction artefact, snake-like pattern on the long-axis view, absence of tendon fibres at the attachment site, peritendinous effusion	65	54.16
Non-retracted (≤ 8 cm)	Type 3a (Fig. 3)	EBA intact or hypertrophied, with minor tendon retraction	10	8.33
Retracted (>8 cm)	Type 3b (Fig. 4)	EBA tear with marked DBBT retraction, no evidence of EBA continuity	55	45.83

Distal biceps: Type 1b pathology




Distal biceps: Type 1b pathology




MRI classification system	Ultrasound classification system	Ultrasound description	Number of cases	Relative (%)
No tear	Type 1	Thickened, hypoechogenic, loss of the fibrillary pattern, tendon continuity to the attachment retained, no evidence of torn fibres on the static and dynamic examinations	20	16.66
	Type 1a (Fig. 7)	It particularly affects one of the two DBBT components (short head or long head)	4	3.33
	Type 1b (Fig. 8)	It affects the full thickness of both DBBT components (short head and long head)	16	13.33
Partial tear	Type 2	DBBT thinning, irregular contour, anechoic appearance, partial tendon discontinuity, peritendinous effusion	35	29.16
	Type 2a (Fig. 9)	Tendon discontinuity at the attachment site, thinning $\leq 50\%$ of the total thickness	18	15.00
	Type 2b (Fig. 5)	Tendon discontinuity at the attachment site, $>50\%$ of the total thickness	14	11.66
High grade (>50 tear)	Type 2c (Fig. 6)	Tendon discontinuity at the attachment site that affects the full thickness of a single DBBT component (short head or long head), which becomes retracted. A refraction artefact at the level of the tendinous stump is observed	3	2.50
	Type 3	Proximal hypertrophic tendon stump with a refraction artefact, snake-like pattern on the long-axis view, absence of tendon fibres at the attachment site, peritendinous effusion	65	54.16
Non-retracted (≤ 8 cm)	Type 3a (Fig. 3)	EBA intact or hypertrophied, with minor tendon retraction	10	8.33
Retracted (>8 cm)	Type 3b (Fig. 4)	EBA tear with marked DBBT retraction, no evidence of EBA continuity	55	45.83

Distal biceps: Type 2a pathology

LAX




SAX




Distal biceps: Type 2a pathology



LAX




SAX



MRI classification system	Ultrasound classification system	Ultrasound description	Number of cases	Relative (%)
No tear	Type 1	Thickened, hypoechogenic, loss of the fibrillary pattern, tendon continuity to the attachment retained, no evidence of torn fibres on the static and dynamic examinations	20	16.66
	Type 1a (Fig. 7)	It particularly affects one of the two DBBT components (short head or long head)	4	3.33
	Type 1b (Fig. 8)	It affects the full thickness of both DBBT components (short head and long head)	16	13.33
Partial tear	Type 2	DBBT thinning, irregular contour, anechoic appearance, partial tendon discontinuity, peritendinous effusion	35	29.16
Low grade (≤ 50 tear)	Type 2a (Fig. 9)	Tendon discontinuity at the attachment site, thinning $\leq 50\%$ of the total thickness	18	15.00
High grade (>50 tear)	Type 2b (Fig. 5)	Tendon discontinuity at the attachment site, $>50\%$ of the total thickness	14	11.66
	Type 2c (Fig. 6)	Tendon discontinuity at the attachment site that affects the full thickness of a single DBBT component (short head or long head), which becomes retracted. A refraction artefact at the level of the tendinous stump is observed	5	2.50
Complete tear	Type 3	Proximal hypertrophic tendon stump with a refraction artefact, snake-like pattern on the long-axis view, absence of tendon fibres at the attachment site, peritendinous effusion	65	54.16
Non-retracted (≤ 8 cm)	Type 3a (Fig. 3)	EBA intact or hypertrophied, with minor tendon retraction	10	8.33
Retracted (>8 cm)	Type 3b (Fig. 4)	EBA tear with marked DBBT retraction, no evidence of EBA continuity	55	45.83






Distal biceps: Type 2b pathology







MRI classification system	Ultrasound classification system	Ultrasound description	Number of cases	Relative (%)
No tear	Type 1	Thickened, hypoechogenic, loss of the fibrillary pattern, tendon continuity to the attachment retained, no evidence of torn fibres on the static and dynamic examinations	20	16.66
	Type 1a (Fig. 7)	It particularly affects one of the two DBBT components (short head or long head)	4	3.33
	Type 1b (Fig. 8)	It affects the full thickness of both DBBT components (short head and long head)	16	13.33
Partial tear	Type 2	DBBT thinning, irregular contour, anechoic appearance, partial tendon discontinuity, peritendinous effusion	35	29.16
Low grade (≤ 50 tear)	Type 2a (Fig. 9)	Tendon discontinuity at the attachment site, thinning $\leq 50\%$ of the total thickness	18	15.00
High grade (>50 tear)	Type 2b (Fig. 5)	Tendon discontinuity at the attachment site, $>50\%$ of the total thickness	14	11.66
	Type 2c (Fig. 6)	Tendon discontinuity at the attachment site that affects the full thickness of a single DBBT component (short head or long head), which becomes retracted. A refraction artefact at the level of the tendinous stump is observed	3	2.50
Complete tear	Type 3	Proximal hypertrophic tendon stump with a refraction artefact, snake-like pattern on the long-axis view, absence of tendon fibres at the attachment site, peritendinous effusion	65	54.16
Non-retracted (≤ 8 cm)	Type 3a (Fig. 3)	EBA intact or hypertrophied, with minor tendon retraction	10	8.33
Retracted (>8 cm)	Type 3b (Fig. 4)	EBA tear with marked DBBT retraction, no evidence of EBA continuity	55	45.83

Distal biceps: Type 2c pathology





Distal biceps: Type 2c pathology






MRI classification system	Ultrasound classification system	Ultrasound description	Number of cases	Relative (%)
No tear	Type 1	Thickened, hypoechogenic, loss of the fibrillary pattern, tendon continuity to the attachment retained, no evidence of torn fibres on the static and dynamic examinations	20	16.66
	Type 1a (Fig. 7)	It particularly affects one of the two DBBT components (short head or long head)	4	3.33
	Type 1b (Fig. 8)	It affects the full thickness of both DBBT components (short head and long head)	16	13.33
Partial tear	Type 2	DBBT thinning, irregular contour, anechoic appearance, partial tendon discontinuity, peritendinous effusion	35	29.16
Low grade (≤ 50 tear)	Type 2a (Fig. 9)	Tendon discontinuity at the attachment site, thinning $\leq 50\%$ of the total thickness	18	15.00
High grade (>50 tear)	Type 2b (Fig. 5)	Tendon discontinuity at the attachment site, $>50\%$ of the total thickness	14	11.66
	Type 2c (Fig. 6)	Tendon discontinuity at the attachment site that affects the full thickness of a single DBBT component (short head or long head), which becomes retracted. A refraction artefact at the level of the tendinous stump is observed	3	2.50
Complete tear	Type 3	Proximal hypertrophic tendon stump with a refraction artefact, snake-like pattern on the long-axis view, absence of tendon fibres at the attachment site, peritendinous effusion	65	54.16
Non-retracted (≤ 8 cm)	Type 3a (Fig. 3)	EBA intact or hypertrophied, with minor tendon retraction	10	8.33
Retracted (>8 cm)	Type 3b (Fig. 4)	EBA tear with marked DBBT retraction, no evidence of EBA continuity	55	45.83



Distal biceps: Type 3a pathology






MRI classification system	Ultrasound classification system	Ultrasound description	Number of cases	Relative (%)
No tear	Type 1	Thickened, hypoechogenic, loss of the fibrillary pattern, tendon continuity to the attachment retained, no evidence of torn fibres on the static and dynamic examinations	20	16.66
	Type 1a (Fig. 7)	It particularly affects one of the two DBBT components (short head or long head)	4	3.33
	Type 1b (Fig. 8)	It affects the full thickness of both DBBT components (short head and long head)	16	13.33
Partial tear	Type 2	DBBT thinning, irregular contour, anechoic appearance, partial tendon discontinuity, peritendinous effusion	35	29.16
Low grade (≤ 50 tear)	Type 2a (Fig. 9)	Tendon discontinuity at the attachment site, thinning $\leq 50\%$ of the total thickness	18	15.00
High grade (>50 tear)	Type 2b (Fig. 5)	Tendon discontinuity at the attachment site, $>50\%$ of the total thickness	14	11.66
	Type 2c (Fig. 6)	Tendon discontinuity at the attachment site that affects the full thickness of a single DBBT component (short head or long head), which becomes retracted. A refraction artefact at the level of the tendinous stump is observed	3	2.50
Complete tear	Type 3	Proximal hypertrophic tendon stump with a refraction artefact, snake-like pattern on the long-axis view, absence of tendon fibres at the attachment site, peritendinous effusion	65	54.16
Non-retracted (≤ 8 cm)	Type 3a (Fig. 3)	EBA intact or hypertrophied, with minor tendon retraction	10	8.33
Retracted (>8 cm)	Type 3b (Fig. 4)	EBA tear with marked DBBT retraction, no evidence of EBA continuity	55	45.83


Distal biceps: Type 3b pathology








Bicipitoradial bursitis




Supination




Pronation

- Surrounds the distal biceps tendon
- Situated between the tendon and the radial tuberosity
- During forearm pronation the bursa is compressed between the radial tuberosity, which becomes posterior, and the distal tendon of the biceps
- Does not communicate with the joint cavity
- May communicate with the interosseous bursa of the elbow

Bicipitoradial bursitis









marc@sonoskills.com *[questions!]*



[Linkedin.com/in/marcsonoskills](https://www.linkedin.com/in/marcsonoskills) *[questions!]*



sonoskills.com



ultrasoundcases.info



SonoSkills WhatsApp +31403041587



[Linkedin.com/company/sonoskills](https://www.linkedin.com/company/sonoskills)



[@sonoskills](https://twitter.com/sonoskills)



[Facebook.com/sonoskills](https://www.facebook.com/sonoskills)



[Instagram.com/sonoskills](https://www.instagram.com/sonoskills)

Live Demonstration



Shelley Guenther, CRGS, CRCS

Sonographer | Clinical Marketing Manager



What additional information would you like?

Interactive Poll

www.clarius.com/msk

www.clarius.com/demo

www.clarius.com/classroom



Clarius L15 HD3

- 30% Smaller & More Affordable
- Wireless Freedom
- High-Definition Imaging
- Easy App for iOS & Android
- Advanced MSK Package
- Clarius Voice Controls
- Clarius Cloud Storage
- Clarius Live Telemedicine
- Unlimited Users

Questions



Marc Schmitz, MSc
Founder & CEO of Sonoskills



Shelley Guenther
Sonographer - @pocus_shelley



Thank you!