

# Clarius Cart

The Clarius Cart gives you both the utility of a cart-based system and the freedom of personal ultrasound. Roll it into a room and start scanning without needing to plug it in.



## Rotating tablet mount

Scan comfortably from any position. The holder fits a compatible iPad\*.



## Ergonomic handle

Easily move the Clarius Cart from room to room.



## Locked storage

Keep equipment secure in a busy hospital. The locked storage container can hold up to 7 Clarius scanners.

## Working space for procedures

An easy-to-clean space to place tools for procedures.




## Charging station

When plugged in, the charging station can charge up to 5 Clarius batteries and 2 smart devices at the same time. The power cable can be tucked away when not in use.



## Wireless ultrasound

Scan using the Clarius Cart without worrying about cords or finding an outlet. Clarius scanners are wireless and connect directly to your phone or tablet.



\* **Compatible iPads:** iPad 2, 3, 4, iPad Air and iPad Air 2. iPad not included.